

Rad-67™ Spot-check Pulse CO-Oximeter®


Featuring a suite of upgradeable rainbow SET™ measurements alongside Masimo SET® Measure-through Motion and Low Perfusion™ pulse oximetry



Oxygen Saturation*



Pulse Rate*



Perfusion Index



Pleth Variability Index



Total Hemoglobin



Methemoglobin**

* Masimo SET® Measure-through Motion technology includes SpO2 and PR

** SpMet is optionally available on Rad-67

More Than a Conventional Pulse Oximeter



Spot-check results displayed with **signal quality indicators**




Label spot-check measurements with **unique patient identifiers** for convenient **historical data review** directly on the device




Measure **SpO2** and **SpHb** with the **rainbow® DCI®-mini sensor** for adult and pediatric patients




Masimo SET® Combined with Next Generation Spot-check SpHb Technology enables:



Measure SpHb, SpO2, pulse rate (PR), and perfusion index (Pi) using the rainbow® DCI-mini sensor **




On-screen guidance to automate workflow



Results displayed in as few as 30 seconds

Features



Intuitive touchscreen allows users to quickly navigate the user interface with finger gestures



Redesigned sensor connector port with a slim profile design provides tactile feedback upon proper connection



HD Display

- Bright LCD, color display
- Automatic low power mode to conserve power

Auto-Brightness

- Ambient light sensor automatically adjusts screen brightness to optimize visibility

Rechargeable Battery

- Li-ion Battery
- Up to 6 hours battery life¹
- 6 hours charging time

Rad-67 Specifications

PHYSICAL CHARACTERISTICS	COMPLIANCE
Weight.....0.37 kg (0.81 lbs)	Safety Standard(s).....ANSI/AAMI ES 60601-1, CAN/CSA C22.2 No. 60601-1, IEC/EN 60601-1, 3rd Ed.
Dimensions.....19.4 cm x 8.2 cm x 2.4 cm (Approx. 7.5" x 3" x 1")	Pulse Oximeter Standard(s).....ISO 80601-2-61
ENVIRONMENTAL	IEC Standard(s).....EN 60601-1-2, Class B
Operating Temperature.....0-35° C (32-95° F)	Type of Protection.....Class II (Internally Powered)
Atmospheric Pressure.....540-1,060 mBar	Degree of Protection.....Type BF, Defib Proof-Applied Part
Operating Humidity.....10-95%, non-condensing	Mode of Operation (per IEC 60601-1).....Continuous Operation
ORDERING INFORMATION	Enclosure Degree of Protection.....IPX4
Rad-67.....PN 9794	

¹This represents approximate runtime at the lowest indicator brightness and wireless functionality powered off using a fully charged battery.

SpHb and SpMet monitoring with Rad-67 are not intended to replace laboratory blood testing. Blood samples should be analyzed by laboratory instruments prior to clinical decision making.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

